

# Elan2

The world's smallest Liquid Nitrogen Generators





"Convenient, reliable, safe, on demand liquid Nitrogen at a touch of a button"

### **Key Advantages:**

- Safe, clean, on-hand Liquid Nitrogen directely in your office/ lab
- + Extracting N2 directly form the air around you
- + Also suitable for HpGe detectors
- Userfriendly interface
- Minimized user intervention
- + Designed with the highest level of user safety in mind
- Improved overall efficiency, uptime and production by integrated self-management system and error checking
- Very clean Liquid Nitrogen —> has not been handled by industrial equipment and transport trucks
- + Ideal solution for research applications, clean room



# MAKE YOUR OWN LIQUID NITROGEN

### **Specifications:**

Voltage: 115 VAC / 230 VAC Generation Rate: approx. 250 ml / h (6—8 L / day)

Phase: Single Phase Storage Capacity: internal 1.0 L Storage

Frequenzy: 50 Hz / 60 Hz Dimensions: Generator: 15" x 15" x 37" (130 lb)

Purity: 98.6 % (< 1 % oxygen) Compressor: 18" x 18" x 15" (60 lb)



### "undemanding systems that deliver the best of results"

### **Operation Environment Requirements:**

- Well ventilated, semi-clean environment
- Room should be free of dust and exhaust
- Temperature should not exceed 77 F
- Humidity: < 80 % rel humidity</li>
- Air Flow: 50 l / min (< 6 ppm moisture)</li>
- Electricity: 2 standard wall outlets
- Power: 900 watts

#### Maintenance:

- An easily accessible air filter is the only required maintenance (approx. every 1—3 months)
- Air filter can be replaced by the customer on-site and with the unit running
- Company MMR-technologies offers a 24 months warranty for the ELAN2 Systems



# MAKE YOUR OWN LIQUID NITROGEN

### **Specifications:**

Voltage: 115 VAC / 230 VAC Generation Rate: approx. 250 ml / h (6—8 L / day)

Phase: Single Phase Storage Capacity: internal 1.0 L Storage

Frequenzy: 50 Hz / 60 Hz + 20.0 L external dewar

Purity: 98.6 % (< 1 % oxygen) Dimensions: Generator: 15" x 15" x 37" (130 lb)

Compressor: 18" x 18" x 15" (60 lb)



# **Denise Fröhlich**

Sales

Cryoandmore Budzylek GbR

Fuggerstraße 9a

41468 Neuss

+ 49 (0) 2131 – 66 46 339

+ 49 (0) 2131 – 66 46 341

denise@cryoandmore.de

# **ALWAYS AT YOUR SERVICE**

